

FAST
TEST



Canine & Feline Giardia Antigen

BF-CFGAG-10
10 TESTS





Result Interpretation



Negative Result

("C")Line: Dark & ("T") Line: Absent
Absence of Giardia in the sample.
Or undetectable levels of Giardia in the sample.



Positive Result

("C")Line: Dark & ("T") Line: Dark
Presence of Giardia in the sample.

Invalid Result

("C")Line: Absent & ("T") Line: Absent
("C")Line: Absent & ("T") Line: Dark
The sample should be re-tested.

