

FAST  
TEST



# Canine

**LSH Ab - ANA Ab**

**EHR Ab - CHW Ag**



BF-CLAEH-10  
10 TESTS

# Result Interpretation



## Negative Result

("C")Line: Dark & ("T") Line: Absent  
Absence of LSH Ab, ANA Ab, EHR Ab, and CHW Ag in the sample.  
Or undetectable levels of LSH Ab, ANA Ab, EHR Ab, and CHW Ag in the sample.



## Positive Result

("C")Line: Dark & ("T") Line: Dark  
Presence of LSH Ab, ANA Ab, EHR Ab, and CHW Ag in the sample.

## Invalid Result

("C")Line: Absent & ("T") Line: Absent  
("C")Line: Absent & ("T") Line: Dark  
The sample should be re-tested.

