

FAST
TEST



Canine CPV & CCV Antigen

BF-CPCAG-10
10 TESTS




Result Interpretation



Negative Result

("C")Line: Dark & ("T") Line: Absent
Absence of CPV & CCV in the sample.
Or undetectable levels of CPV & CCV in the sample.



Positive Result

("C")Line: Dark & ("T") Line: Dark
Presence of CPV & CCV in the sample.



Invalid Result

("C")Line: Absent & ("T") Line: Absent
("C")Line: Absent & ("T") Line: Dark
The sample should be re-tested.

