

FAST
TEST



Canine CPV-CCV-Giardia Antigen

BF-CPCGA-05
5 TESTS




Result Interpretation




Negative Result

(“C”)Line: Dark & (“T”) Line: Absent
Absence of CPV, CCV & Giardia in the sample.
Or undetectable levels of CPV, CCV & Giardia in the sample.



Positive Result

(“C”)Line: Dark & (“T”) Line: Dark
Presence of CPV, CCV & Giardia in the sample.



Invalid Result

(“C”)Line: Absent & (“T”) Line: Absent
(“C”)Line: Absent & (“T”) Line: Dark
The sample should be re-tested.

