

FAST
TEST



E. canis

Antibody

BF-ECAB-10
10 TESTS




Result Interpretation




Negative Result

("C")Line: Dark & ("T") Line: Absent
Absence of *E. canis* antibodies in the sample.
Or undetectable levels of *E. canis* antibodies in the sample.



Positive Result

("C")Line: Dark & ("T") Line: Dark
Presence of *E. canis* antibodies in the sample.



Invalid Result

("C")Line: Absent & ("T") Line: Absent
("C")Line: Absent & ("T") Line: Dark
The sample should be re-tested.

