

FAST
TEST



Feline FIV Antibody

BF-FIAB-10
10 TESTS





Result Interpretation



Negative Result

("C")Line: Dark & ("T") Line: Absent
Absence of FIV antibodies in the sample.
Or indetectable levels of FIV antibodies in the sample.



Positive Result

("C")Line: Dark & ("T") Line: Dark
Presence of FIV antibodies in the sample.

Invalid Result

("C")Line: Absent & ("T") Line: Absent
("C")Line: Absent & ("T") Line: Dark
The sample should be re-tested.

