

FAST
TEST



Feline

FIV Antibody

FeLV Antigen

BF-FILAT-10
10 TESTS




Result Interpretation




Negative Result

("C")Line: Dark & ("T") Line: Absent
Absence of FIVAb & FeLV Ag in the sample.
Or undetectable levels of FIV Ab & FeLV Ag in the sample.



Positive Result

("C")Line: Dark & ("T") Line: Dark
Presence of FIVAb & FeLV Ag in the sample.



Invalid Result

("C")Line: Absent & ("T") Line: Absent
("C")Line: Absent & ("T") Line: Dark
The sample should be re-tested.

