

FAST
TEST



Feline Toxoplasma Antibody


BF-FTAB-10
10 TESTS

Result Interpretation



Negative Result

("C")Line: Dark & ("T") Line: Absent
Absence of toxo antibodies in the sample.
Or undetectable levels of toxo antibodies in the sample.



Positive Result

("C")Line: Dark & ("T") Line: Dark
Presence of toxo antibodies in the sample.

Invalid Result

("C")Line: Absent & ("T") Line: Absent
("C")Line: Absent & ("T") Line: Dark
The sample should be re-tested.

