

FAST  
TEST



# AIV Antigen

BF-AAG-10  
30 TESTS




# Result Interpretation




## Negative Result

("C")Line: Dark & ("T") Line: Absent  
Absence of AIV in the sample.  
Or undetectable levels of AIV charge in the sample.



## Positive Result

("C")Line: Dark & ("T") Line: Dark  
Presence of AIV in the sample.



## Invalid Result

("C")Line: Absent & ("T") Line: Absent  
("C")Line: Absent & ("T") Line: Dark  
The sample should be re-tested.

