

FAST
TEST



Canine Total IgE



BF-CTIGE-10
10 TESTS

Result Interpretation



Negative Result

("C")Line: Dark & ("T") Line: Absent
Absence of IgE antibodies in the sample.
Or undetectable levels of IgE antibodies in the sample.



Positive Result

("C")Line: Dark & ("T") Line: Dark
Presence of IgE antibodies in the sample.

Invalid Result

("C")Line: Absent & ("T") Line: Absent
("C")Line: Absent & ("T") Line: Dark
The sample should be re-tested.

