

FAST
TEST



Feline FCoV Antigen


BF-FCAG-10
10 TESTS

Result Interpretation




Negative Result

("C")Line: Dark & ("T") Line: Absent
Absence of FCoV in the sample.
Or undetectable levels of FCoV in the sample.



Positive Result

("C")Line: Dark & ("T") Line: Dark
Presence of FCoV in the sample.



Invalid Result

("C")Line: Absent & ("T") Line: Absent
("C")Line: Absent & ("T") Line: Dark
The sample should be re-tested.

