

FAST
TEST



Feline FCoV Ag & Ab

BF-FCAGB-05
5 TESTS



Result Interpretation




Negative Result

("C")Line: Dark & ("T") Line: Absent

Absence of FCoV antigens & antibodies in the sample.


Or undetectable levels of FCoV antigens & antibodies in the sample.



Positive Result

("C")Line: Dark & ("T") Line: Dark

Presence of FCoV antigens & antibodies in the sample.



Invalid Result

("C")Line: Absent & ("T") Line: Absent

("C")Line: Absent & ("T") Line: Dark

The sample should be re-tested.

